
Student name:

Code number:

GENETICS

Lab 2

(2014 - 2015)

Please read and make sure you understand the following instructions and knowledge before you go on.

Revised from Lecture 2:

- The Cell cycle
- Mitosis phases

1. CELL CYCLE AND MITOSIS

Objective:


Upon completion of this lab, the students should be able to:

1. Understand the cell cycle process and different stages of mitosis.
2. Know the importance of feulgen technique in cytological research.
3. Prepare his/her own root tip squashing technique to observe mitotic divisions in meristematic cells of Faba bean roots' tips.
4. Evaluate the quality of the stain in relation to hydrolysis.
5. Visualize different phases of mitosis on broad beans (*Vicia faba*) root tip.
6. Identify the various stages of mitosis in the root cells.
7. Sketch each mitotic stage and summarize the main events of each stage.
8. Describe the events during each phase of mitosis (number and structure).
9. Value the mitotic index and the duration of each stages in the cell cycle.

Eukaryotic organisms carry out mitosis throughout their entire life to grow, develop, and asexually reproduce (some cases), as well as to renew the old and damaged cells.

CELL CYCLE

The cell cycle is the series of events that take place in a cell leading to its division and duplication (replication). In cells without a nucleus (prokaryotic), the cell cycle occurs via a process termed binary fission. In cells with a nucleus (eukaryotes), the cell cycle can be divided in three periods: **interphase**—during which the cell grows, accumulating nutrients needed for mitosis and duplicating its DNA—and **the mitotic (M) phase**, during which the nucleus of the cell splits itself into two, forming two "daughter cells" that is completely divided in the final phase, **cytokinesis**.



Interphase

The cell cycle occurs from the completion of one division until the completion of the next division. It involves 2 phases: Interphase (G₁, S and G₂) and Mitosis (M) followed by Cytokinesis (C). The period between M and S is called G₁ stage and that between S and M is G₂ stage. The cell spends 90% of its time in Interphase and only 10% in Mitosis but, the duration of each phase and stage in eukaryotic cells depends on the cell type.

Mitosis

Cell growth stops at this stage and cellular energy is focused on the orderly division into two daughter cells.


The relatively brief *M phase* consists of nuclear division (karyokinesis). It is relatively short period of cell cycle. M phase is complex and highly regulated. The sequence of events is divided into phases, corresponding to the completion of one set of activities and the start of the next. These phases are sequentially known as prophase, metaphase, anaphase and telophase.

Cytokinesis


It is not part of mitosis but is an event that directly follows mitosis in which cytoplasm is divided into two daughter cells, each with the same number and same quality of chromosomes, but with different structures

(dyads to monads). The cell then enters interphase - the interval between mitotic divisions.

Mitosis - *Allium* Root Tip




Interpahase




Prophase




Metaphase



Anaphase



Telophase




Late Telophase
and Cytokinesis

FEULGEN SQUASH TECHNIQUE FOR MITOSIS

Feulgen stain is undoubtedly the most widely used staining technique discovered by Robert Feulgen. It is a semi-quantitative technique used in cyto-histology for investigating chromosomes and its DNA (DNA-specific reaction) in dividing cells depending on acid hydrolysis of DNA and be visualized with optical microscopy.

In plants, mitotic cell division mainly takes place in special regions called **meristems**. They are either present in Shoot apex or axillary buds or root tips of the plants for development and growth.



MATERIALS AND METHOD

Material

1. Previously prepared broad bean (*Vicia faba*) root tips.
2. Light microscope with 10x and 40x objectives.
3. Water bath
4. Lab coat
5. Slides and coverslips
6. Dry, clean test tubes in rack.
7. Paint brush
8. Pencil eraser
9. Razor blade
10. Paper towels

Chemicals preparation

- 1. 1N Hydrochloric acid (HCl):** If your stock concentration came in 12M (smocking), then using the following equation:

$$N V = N_1 V_1$$

$$\text{So, } V = (N_1 V_1) / N$$

Where,

V = volume of your stock concentration, N = concentration of your stock

V₁ = volume of your final preparation, N₁ = concentration of your final


$$V = (1 \text{ N} \times 100 \text{ ml}) / (12 \text{ N}) = 8.3 \text{ ml}$$

- 2. 45% acetic acid:** 45 ml of glacial acetic acid rises up to 100 ml with distilled water.
- 3. Feulgen stain (Schiff's reagent):** one gram basic fuchsin was dissolved in 200 ml boiled distilled water, shake well and cool to 50-55°C. Add 30 ml 1N HCl to the stain then 3 g potassium metabisulphide (K₂S₂O₅). Allow the stain to bleach for 24 hr in a dark stoppered bottle. Shake well, filter with charcoal and stored the bottles in dark.

Procedure


1. Dry broad beans (uniform in size and shape) were soaked in tap water for 24 hr. soaking water must be changed every 12 hr.

Why? to avoid anoxia of cells on divisions.



2. Seeds were sown vertically in irrigated loose sand or sawdust for 3 days. **Why loose sand or sawdust was used?** So, the root will not rupture by harvesting.

Germination of bean seed



3. Healthy straight roots (about 2 cm) were collected and washed with distilled water to remove excess sand or sawdust. **Why we take the root tip?** To see mitotic divisions as it consist of meristematic dividing cells.



4. The root tips (0.5-1 cm) were cut off and transfer immediately the root tips in Fixative solution (acetic acid : alcohol =1: 3) for 24 hr, then store them in 70% ethanol in refrigerator. **Why we do this process?** The combination of these two chemicals leaves the chromosomes in a highly precipitated form suitable for staining and microscopic study by preserving their structural and/or chemical components i.e. stop the cell division and growth.



5. Wash the needed roots in test tubes 5 times for 5 min. **Why?** to remove excess fixative.
6. Separate the tips into 3 different labeled tubes.
7. Hydrolyze the root tips in 1N HCl (test tubes in water bath) at constant 60°C for 7, 8 and 9 min in separate tubes. **Why the**

tissues were hydrolyzed? 1.To separate the cells from each other to see the division clearly. 2.To break the purine-base glycoside linkage yielding free aldehyde group from deoxypentose sugar along DNA.



8. Pour off the fixative and wash the hydrolyzed roots in distilled water 3 times to get rid of the acid (check that the odor of acid is removed).
9. Using the paint brush to transfer the tissue to tubes filled with 2ml feulgen reagent (Schiff reagent). **Be careful that the roots are now very fragile.**
10. Stain for $\frac{1}{2}$ hour in the lab with leucobasic fuchsin in dry clean tubes. **What happened?**

The Schiff reagent is prepared by adding 1N HCl (acid) and potassium meta-bisulphite (salt) to an aqueous solution of basic fuchsin. This addition causes the liberation of free SO_3 which reduces the purple color of the basic fuchsin into a colorless solution (leucobasic fuchsin). The reagent reacts with the free aldehydic groups (CHO) from the deoxypentose sugar of DNA producing a DNA specific insoluble magenta colour compound.

Note: Take care this dye will stain hands and clothes.



11. Remove the root tip from the stain and transfer it to a clean slide.



12. Cut the dark purple color of the **meristematic region** with sharp razor blade.





13. Squash the meristematic region of the root tip with the pencil eraser and add 1-2 drops of 45% acetic acid onto the slide. **Why using acetic acid on squashing?** To easily spread the cells in a very thin layer to see the divisions clearly.

14. Place a coverslips over the tissue.

Note: Try to avoid air bubbles under the coverslip to prevent drying the tissue.

15. Press down carefully onto the coverslip with pencil eraser to spread the cells in a very thin layer to see the divisions. The tapping should be gentle as the hard way will split the cover.


Note: The coverslips should not slip or break.



16. Use the paper towels or blotting paper to remove excess acetic acid over the slide and cover slip.

17. Examine under low power (10x objective) to found areas with cells undergoing mitosis.

18. Examine under high power (40 x objective) and find all the stages.




19. After you have set up your own slides, your instructor will adjust your best slides on the main lab bench pointing on the stages of mitosis.


Observation:

1. The material was spread into one layer of flattened rectangular cells.

Note: If not you will not be able to see divisions. You will see only overlapped cells.



2. All the Cell cycle stages were present: Interphase, mitosis and cytokinesis.



- a. Metaphase, b. Anaphase, c. Telophase, d. Prophase and e. Interphase
3. All the cells appear rectangular.

Note: If not (the cell in spindle form) your coverslip must have moved.


4. The best cell stain must have **purple color**.

Why?

Optimum hydrolysis liberates adequate purine and aldehydic groups from the deoxypentose sugar of DNA to react with the staining reagent giving rise to the maximum results of the dye (**purple color**).

But if **faint color**, it may be the effect of either:

- less hydrolysis due to the few liberation of purine and aldehydic groups
- **or** excess hydrolysis due to destructive effect on the structure of nucleic acid by removing both histones (protein associated with DNA) and apurinic acids (deoxyribose without the base) .



WORK SHEET

1. Examine your own slide under microscope. Scan your slide to find all the mitotic stages and draw them below. Indicate the chromosome number and structure of each stage.

	Interphase	Mitosis	
		Prophase	Metaphase
Draw			
chromosome number			
chromosome structure			
Draw		Anaphase	Telophase
chromosome number			
chromosome structure			

2. Record the total number of cells that you can see in the field of your microscope. Calculate the percentage of each stage and its duration in the cell cycle. Where the cycle duration=1440 min (divide every 24 hr).

The duration of Stage in the cell cycle (minutes) = (number of cell in stage / total number of cells) x 1440

Stage of cell cycle	# cells in the stage	Duration of the stage (min)
Interphase		
Prophase		
Metaphase		
Anaphase		
Telophase		
Total cells cells

Note: The duration of each phase and stage in eukaryotic cells depends on the cell type.

3. Quantify the cell division of the broad bean in your slide by using the mitotic index:

Mitotic index = number of cell in mitosis / total number of cells


.....

Note: The **mitotic index** is simply a measurement to determine the percentage of cells undergoing mitosis. Mitosis is the division of somatic cells when genetic information from one single cell is equally dispersed

into two daughter cells. Durations of the cell cycle and mitosis vary in different cell types. An elevated mitotic index indicates more cells are dividing; and thus obvious in cancer cells. The mitotic index may be elevated during necessary processes to life, such as the normal growth of plants or animals, as well as cellular repair the sire of an injury.

Student name:**Code number:**

STUDENT'S ASSINMENT
Give it to your laboratory instructor



Q1. Why RNA is not hydrolyzed by the HCl treatment?

.....

.....


.....

.....

.....

.....

Q2. Examine the provided photo of the L.S. in root apical meristem. Are most of the cells dividing or is the majority in interphase?



.....

.....

.....

.....

.....

Q3. Data were collected from a microphotograph of onion and bean root tips when they are dividing, showing the following data:

Stage of cell	# cells in the stage
---------------	----------------------

cycle	Onion	Bean
Interphase	64	3
Prophase	18	18
Metaphase	3	64
Anaphase	6	9
Telophase	9	6
Total number of cells	100 cells	100 cells

- Calculate the percentage of each stage and its duration in the cell cycle as min for both plants. Where the cycle duration=1440 min (divide every 24 hr).
- Calculate the mitotic index of both plants and comment on your data.

This image shows a full page of white paper with horizontal dotted lines, typical of primary-ruled notebook paper. The lines are evenly spaced and run across the entire width of the page. There are no margins, text, or other markings present.

[illegible]

Q4. Choose the correct answer:

1. The stage of the cell cycle where each chromosome is composed of two chromatids in preparation for mitosis.

- a. G1 b. S c. M d. G2

2. During which stage of mitosis do the centromeres split?

- a. Prophase b. Metaphase c. Anaphase d. Telophase

3. During which stage of mitosis does the nuclear envelope begin to disappear?

- a. Prophase b. Metaphase c. Anaphase d. Telophase

4. The pinching off of the cell membrane that creates two new cells (after mitosis) is called

- a. Interphase b. Mitosis c. Cytokinesis d. Meiosis

5. A cell with 10 chromosomes undergoes mitosis. How many daughter cells are created? ____ Each daughter cell has ____ chromosomes.

- a. 2, 10 b. 10, 2 c. 1, 10 d. 2, 20

6. If a diploid cell has the genotype Pp, what will be the genotype of the two daughter cells after mitosis?

- a. PP and pp b. P and p c. Both PP d. Both pp e. Both Pp