

Determination of HCl sample




1- Principle

HCl

Directly determined, by titration with 0.1 N NaOH using any indicator ?

End point: Ph.Ph. colorless → **Pink** (PH range 8–10)


M.O **Red** → **Orange** (PH range 3.3.–4.2)



pH= 7
Neutral

2- Procedure

In a conical flask



Titrant
(filled in
burette)

0.1 N NaOH


End Point: First
persistent **pink**

Record reading of the
burette


Indicator
+
Sample

Add 10 drops Ph.Ph.

Add 10 ml HCl sample (by bulb pipette)



In a conical flask



**Titrant
(filled in
burette)**

0.1 N NaOH


End Point: Orange

**Record reading of the
burette**

**Indicator
+
Sample**

Add 2 drops M.O.

Add 10 ml HCl sample (by bulb pipette)



3. Calculation

1ml standard \equiv ?? g sample

1mL 0.1 N NaOH \equiv F g HCl

Equivalent Factor

Balanced equation $1 \text{ NaOH} + 1 \text{ HCl} \rightarrow \text{NaCl} + \text{H}_2\text{O}$

Std (left) \equiv Sample(right)

1 mole NaOH \equiv 1 mole HCl

1L 1M NaOH \equiv Mwt. HCl

1L $\frac{1M}{1}$ NaOH \equiv $\frac{\text{Mwt.}}{1}$ HCl

1L 1N NaOH \equiv $\frac{\text{Mwt.}}{1}$ HCl

1L 0.1N NaOH \equiv $\frac{\text{Mwt.}}{1 \times 10}$ HCl

1ml 0.1N NaOH \equiv $\frac{\text{Mwt.}}{1 \times 10 \times 1000}$ HCl

$1M = \frac{1 \text{ mole}}{1L} \rightarrow 1 \text{ mole} = 1L \text{ } 1M$

$N = \frac{\text{Equivalent weight}}{L}$

$N = \frac{M}{\text{no. of OH}^-}$

Divide both sides by 10

Divide both sides by 1000

$F = \frac{36.5}{1 \times 10 \times 1000} = 0.00365 \text{ g}$

$$\text{Conc. of HCl} = \frac{\text{mls} \times f \times F \times 1000}{10} = \text{g/L}$$



Determination of Acetic acid sample




1 – Principle

- ▶ CH_3COOH is Weak acid which is directly titrated against NaOH

weak acid \neq strong base

- ▶ At the end point, the following are present




PH is Slightly Alkaline

ph.ph. indicator

2- Procedure

In a conical flask



Titrant
(filled in
burette)

0.1 N NaOH


End Point: First
persistent **pink**

Record reading of the
burette

Indicator
+
Sample

Add 10 drops Ph.Ph.

Add 10 ml Acetic acid sample
(by bulb pipette)



3- Calculation

1ml standard \equiv ?? g sample



Std (left) \equiv Sample(right)

1 mole NaOH \equiv 1 mole CH₃COOH

1M = 1 mole/1L

1L 1M NaOH \equiv Mwt. CH₃COOH

1L $\frac{1M}{1}$ NaOH \equiv Mwt. CH₃COOH

1L 1N NaOH $\equiv \frac{\text{Mwt. CH}_3\text{COOH}}{1}$

N = Equivalent weight/L

Divide both sides by 10

1L 0.1N NaOH $\equiv \frac{\text{Mwt. CH}_3\text{COOH}}{1 \times 10}$

Divide both sides by 1000

1ml 0.1N NaOH $\equiv \frac{\text{Mwt. CH}_3\text{COOH}}{1 \times 10 \times 1000}$

$$\mathbf{F} = \frac{60}{1 \times 10 \times 1000} = 0.006 \text{ g}$$

$$\text{Conc. of acetic acid} = \frac{\text{mls} \times f \times F \times 1000}{10} = \text{g/L}$$



Thank You

