


# Determination of Sodium Carbonate sample




# Principle

$\text{Na}_2\text{CO}_3$  is salt of strong base & weak acid so can be titrated against strong acid using ph.ph or M.O.



	$\text{CO}_3^{--}$	$\xrightarrow{\text{H}^+}$	$\text{HCO}_3^-$	$\xrightarrow{\text{H}^+}$	$\text{CO}_2 + \text{H}_2\text{O}$
<b>pH</b>	11.7		8.3		3.8
<b>Ph.Ph.</b> (8-10)	Pink		Colorless $\frac{1}{2}$ neutralization		Colorless
<b>M.O.</b> (3.3-4.4)	Yellow		Yellow		Orange

**a) Using M.O. as indicator:**





*Procedure*

- 1) 10 mL sample in a conical flask + 2 dps M.O.
- 2) Titrate against 0.1 N HCl
- 3) End point: **yellow** to **orange**

*Calculation*



$$\text{Concn. of Na}_2\text{CO}_3 = [\text{mls} \times f \times F \times 1000] / 10$$



$$F = 0.0053 \text{ g}$$

## b) Using Ph.Ph. as indicator:



2 Steps





## How to avoid the escape of CO<sub>2</sub>?


1	<b>Dilution</b>	→ To increase solubility of CO <sub>2</sub>
2	<b>Cooling</b>	
3	<b>Slow titration</b>	→ To avoid sudden formation of large amounts of CO <sub>2</sub>
4	<b>Stirring</b>	→ To distribute CO <sub>2</sub> in the solution to allow all CO <sub>2</sub> to react
5	<b>Dipping the nozzle of the burette in the solution surface</b>	



# Procedure

In Beaker (250 ml)

**Glass rod: Once it's put, NEVER removed until the end of titration**



- 1** 10 ml sample + 10 dps Ph.Ph.
- 2** Dilute till 2/3 beaker with distilled water
- 3** Cool very well in ice bath
- 4** Titrate against 0.1N HCl with **dipping nozzle** of burette in the solution with stirring till disappearance of pink color

# Calculation

$$\text{Concn. of Na}_2\text{CO}_3 = [\text{mls} \times f \times F \times 1000] / 10$$



$$1 \text{ mole HCl} = 1 \text{ mole Na}_2\text{CO}_3$$

$$1 \text{ L } 1 \text{ M HCl} = \text{Mwt Na}_2\text{CO}_3$$

$$1 \text{ mL } 0.1 \text{ N HCl} = \text{Mwt Na}_2\text{CO}_3 / 10 \times 1000$$

$$F = 0.0106 \text{ g}$$

# **Determination of Calcium Carbonate sample**





# Uses

- ❑ Source of Calcium in tablets
- ❑ Antacid
- ❑ Chalk



# Types of Titration

**Direct Titration**

**Acid  $\neq$  Base**

**Indirect Titration**

**Sample + Reagent**

**Back Titration**

**Acid  $\neq$  Base**

- 1) Insoluble substances, e.g. ZnO CaCO<sub>3</sub>
- 2) Volatile substances, e.g. formic acid or ammonia.
- 3) Reactions that need heat.
- 4) Slow reactions.

# 1- Principle

## Back Titration

.. Why ?

Because  $\text{CaCO}_3$  is water insoluble

.. How ?

$\text{CaCO}_3$  react with **Known Excess Standard HCl** →  
produce  $\text{H}_2\text{CO}_3$  ( $\text{CO}_2$  and  $\text{H}_2\text{O}$ )


**Excess Unreacted (Remaining) HCl** is **Back Titrated**  
with standard NaOH using Ph.Ph. as indicator till first pink  
color

## N.B.1

### Back titration is for ...

- 1) Insoluble substances, e.g., ZnO or CaCO<sub>3</sub>
- 2) Volatile substances, e.g., formic acid or ammonia.
- 3) Reactions that need heat.
- 4) Slow reactions.

## N.B. 2



## 2- Procedure

Weigh about 0.3g  $\text{CaCO}_3$  powder on glazed paper

Transfer to conical flask

Add 25ml **0.5N** HCl

(with bulb pipette  
or burette)

Boil for 2 min

Why

- Help  $\text{CaCO}_3$  to react with HCl
- To evolve  $\text{CO}_2$

Cool

Why

Ph.Ph. is prepared  
in alcohol

Add 10 drops Ph.Ph.

Titrate against **0.5N** NaOH


End Point: first pink

# 3- Calculation

F


We have 2 standards but **F** is a relation between sample and standards reacting with each other  
i.e. **first standard** that meet the **sample**

HCl



CaCO<sub>3</sub>

Balanced equation



$$\text{Concn. of CaCO}_3 = \frac{(25 \times f_1 - \text{mls} \times f_2) \times F \times 100}{\text{weight of sample}(g)} = \boxed{g\%}$$

$f_1$ : Standardization factor for HCl standard

$f_2$ : Standardization factor for NaOH standard



**Remember:**  
The sample is Solid

2 mole HCl = 1 mole CaCO<sub>3</sub>

1 mL 0.5 N HCl = M.Wt. CaCO<sub>3</sub> / (2 × 1000 × 2)

F = 100 / (2 × 1000 × 2) = 0.025 g



**Thank You**